

INTRODUCING NEW DYNA GRIP™

Powder-Free Latex Exam Grade Gloves

- Silky Smooth Inside for Easy Slip-On
- Specially Textured Surface Creates Excellent Grip, Wet or Dry
- Exclusive Dyna-Flex Latex Formula for Superior Dexterity and Comfort
- Extra Thick, 8 Mil. Powder-Free Latex Construction




GET GRIP!



INTRODUCING NEW DYNA GRIP™

Powder-Free Latex Exam Grade Gloves

Quality Hand Protection From:



Textured Surface for Incredible Gripping Power



Extreme Flexibility



Easy Slip-On Comfort



SPECIFICATIONS

Material: Latex, Powder-Free
Thickness: 8 Mil
Glove Length: 9" With Beaded Cuff
Protein Level: Less than 50 µg/g
Packaging: 100 Gloves By Weight Per Box

Size:	Part No.
Extra Small	650-1000
Small	650-1001
Medium	650-1002
Large	650-1003
Extra Large	650-1004
XX Large	650-1005

FINISHED PRODUCT SPECIFICATION AND STANDARDS:

POWDER FREE LATEX EXAMINATION GLOVES WITH LOW PROTEIN			
510(K) Number: K022088			
Type: Latex 9" 8mil	Color: Natural		
Protein Level: less than 50 µg/g			
Typical characteristic	Overall textured surface/beaded cuff/Ambidextrous	Lubrication	Chlorinated

Size	SMALL	MEDIUM	LARGE	X-LARGE
Length (mm)	240-260	240-260	240-260	240-260
Width (mm)	70-90	85-105	100-120	110-130

STERILE

Size	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
Cuff (mm)	0.12	0.12	0.12	0.12	0.12	0.12
Palm (mm)	0.19	0.19	0.19	0.19	0.19	0.19
Finger (mm)	0.20	0.20	0.20	0.20	0.20	0.20

Size	SMALL	MEDIUM	LARGE	X-LARGE
Average Weight (gm)	8-9	9-10	10-11	11-12

PACKAGING

Item type	Inner box	Carton box	1_20_Container	1_40_Container	Note
Type	100 pieces	1000 pcs	800,000 pcs	1600,000 pcs	

STANDARD AND REQUIREMENTS

Item type		AQL	Level
Testing Methods	ASTM D3578-01		
Standards	ISO2859		
Classifications	Water Leak Test	AQL2.5	G-I
	Major Defects	AQL2.5	G-I
	Minor Defects	AQL4.0	G-I

PHYSICAL PROPERTIES

Item type		Before Aging (min)	After Aging (min)
Tensile Strength	FDA minimum requirement	18MPa	14MPa
	Our Average	18MPa	14MPa
Ultimate Elongation	FDA minimum requirement	650%	500%
	Our Average	650%	500%
Shelf Life		3 years	

Shipper Carton



Dispenser Box



Signal Hill, CA 90755
 1-800-262-0200 • sassafety.com

P/N DynaGripGF001