



ChemDefender® disposable gloves are engineered with a high-grade synthetic rubber chloroprene material that allows for superior durability while providing total flexibility and comfort. Chloroprene, compared to Nitrile, exceeds on strength and chemical resistance. ChemDefender® gloves are made with an extended cuff for added protection and textured fingertips for better grip action.

## 5 mil CHLOROPRENE

- PF**  
POWDER-FREE
- EX**  
EXAM GRADE
- LATEX**  
LATEX-FREE
- TEXTURED GRIP**
- SINGLE USE ONLY**
- EXTENDED CUFF**

### Powder-Free Chloroprene Disposable Gloves

- Powder-free, exam grade chloroprene
- Thickness: 5 mil
- Non-latex
- Superior dexterity and durability
- 12" length for added protection
- Textured finger tips for better grip
- Single use only

Box of 50 Item No.	Size	Case Qty	UPC	ITF10
66591	Small	10 Boxes	781311665910	10781311665917
66592	Medium	10 Boxes	781311665927	10781311665924
66593	Large	10 Boxes	781311665934	10781311665931
66594	X-Large	10 Boxes	781311665941	10781311665948
66595	2X-Large	10 Boxes	781311665958	10781311665955



### Supreme Chemical Resistance\*



\*vs Nitrile after a 15 minute chemical bath



- Competitive Nitrile
- Acetone
- .5 lbs pressure
- 55 mm



- Chem Defender
- Acetone
- 2 lbs pressure
- 141 mm



- Competitive Nitrile
- MEK
- .25 lbs pressure
- 35 mm



- Chem Defender
- MEK
- 1 lb pressure
- 112 mm



**SCAN ME FOR MORE INFO**

#### Additional Specs

Palm Width	
Small	86 ± 4 mm
Medium	95 ± 4 mm
Large	106 ± 4 mm
X-Large	115 ± 4 mm
2X-Large	120 ± 4 mm

#### Tensile Strength

Before Aging (MPa)	min. 14
After Aging (MPa)	min. 14

#### Elongation

Before Aging (%)	500
After Aging (%)	400

#### Length / Thickness (single wall)

Length	290 ± 5 mm
Cuff Thickness	0.09 ± 0.02 mm
Palm Thickness	0.12 ± 0.02 mm
Finger Thickness	0.16 ± 0.02 mm