



POWDERED OR
POWDER-FREE



EXAM GRADE



LATEX-FREE



TEXTURED
GRIP



SINGLE USE
ONLY

Nitrile Disposable Gloves

- Powdered nitrile or powder-free exam grade nitrile
- Thickness: 5 mil
- Non-latex
- Extraordinary strength and puncture resistance
- Lightly textured
- Single use only



Powder-Free

Box of 100 Item No.

6606-20

6607-20

6608-20

6609-20

Size

Small

Medium

Large

X-Large

Powdered

Box of 100 Item No.

6606

6607

6608

6609

6610

Size

Small

Medium

Large

X-Large

2X-Large



Palm Width (mm)

Size	Derma-Lite®	ASTM Standard
Small	85 ± 5	80 ± 10
Medium	95 ± 5	95 ± 10
Large	105 ± 5	111 ± 10
X-Large	115 ± 5	120 ± 10
2X-Large	120 ± 5	120 ± 10

Length / Thickness (single wall)

Dimension	Derma-Lite®	ASTM Standard
Length (mm)	245 ± 10	120 ± 10
Cuff Thickness (mm)	0.09 ± 0.02	N/A
Palm Thickness (mm)	0.10 ± 0.02	0.05
Finger Thickness (mm)	0.16 ± 0.02	0.05

Tensile Strength

	Derma-Lite®	ASTM Standard
Before Aging (MPa)	min. 14	min. 12.5
After Aging (MPa)	min. 14	min. 12.5

Elongation

	Derma-Lite®	ASTM Standard
Before Aging (%)	500	500
After Aging (%)	400	400