



RAVEN® EXAM GRADE NITRILE DISPOSABLE GLOVES

Raven® powder-free exam-grade nitrile disposable gloves are made from advanced proprietary materials, providing excellent performance along with chemical and puncture resistance. With a thickness of 7 mils, these gloves offer premium hand protection and are textured for improved grip, ensuring they can handle even the toughest jobs with confidence. Ideal for various applications including industrial work, automotive tasks, janitorial services, food handling, staining, painting, and working with oil-based paints. The Raven® is the best-selling black nitrile disposable glove in the industry!



INDUSTRIES & APPLICATIONS

- Automotive
- Chemical
- Correctional Facilities
- Dental
- Electronics
- Food Industries
- Healthcare
- Janitorial
- Machinery and Equipment
- Mining
- Oil and Gas
- Plumbing
- Recycling and Waste Disposal
- Tattoo Studios
- Transportation
- Utilities
- Warehousing
- Painting

PRODUCT FEATURES

- Powder-free exam grade nitrile
- Thickness: 7 mil
- Non-latex
- Outstanding strength, wear and dexterity
- Textured for enhanced gripping power
- Passes ASTM D6978 to handle Fentanyl substances
- Made from component materials which comply with FDA regulations for food contact. 21 CFR 170-199





STANDARDS & PERFORMANCE

FDA Regulations	21 CFR 170-199 Food contact
Inspection Level	AQL 10.0 Level S3
Standards	ASTM D6978, ASTM D6319
Storage	Keep out of direct sunlight; store in a cool and dry place.

TECHNICAL DATA

Material	Nitrile
Thickness	7 mil
Length	9.5"
Grade	Exam
Powder	Powder-Free
Color	Black
Texture	Fully
Sizes	S - 2XL
Usage	Single Use
Country of Origin	Malaysia

PHYSICAL PROPERTIES

Length (mm / inches)	245 / 9.6
Palm Thickness (mm/mil)	0.14 / 5.51
Finger Thickness (mm/mil)	0.18 / 7.09
Tensile Strength	14 MPa (min)
Tensile Strength (Aged)	14 MPa (min)
Elongation	500%
Elongation (Aged)	400%

ITEM INFORMATION

Box of 100

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF
66516	Small	100 Gloves per box	9.25 x 5.25 x 3.75	1.48	10 Boxes per case	15.5 x 11.12 x 10.25	16.24	781311665163	10781311665160
66517	Medium			1.58			17.74	781311665170	10781311665177
66518	Large			1.65			18.81	781311665187	10781311665184
66519	X-Large			1.82			19.56	781311665194	10781311665191
66520	2X-Large			1.98			20.35	781311665200	10781311665207

Box of 50

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF
66517-01	Medium	50 Gloves per box	6.25 x 5.75 x 3	0.81	12 Boxes per case	18.75 x 10 x 7	10.76	781311665071	10781311665078
66518-01	Large			0.875			11.13	781311665088	10781311665085
66519-01	X-Large			0.91			12.34	781311665095	10781311665092

Pack of 6

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Inner Qty	Inner Dims (in)	Inner Wt (lb)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF10	ITF20
66511	Large	6 Gloves per pack	9.25 x 5.25 x 0.38	0.105	25 Packs per inner	10.38 x 5.25 x 7.5	2.63	50 Packs per case	10.24 x 9.45 x 5.91	6.04	781311665118	10781311665115	20781311665112
66512	X-Large			0.106			2.65			6.09	781311665125	10781311665122	20781311665129

Read and Understand Instructions and Warnings before using this product.

Raven Safety makes no recommendations as to the use of this product. All safety equipment should be selected in conjunction with an OSHA approved safety program that evaluates your specific hazards. Improper selection and/or use of safety products may cause serious injury.

