



# THICKSTER®



## THICKSTER EXAM GRADE LATEX DISPOSABLE GLOVES

Thickster® gloves offer excellent protection and performance. 14 mils thick, more than double the thickness of standard latex gloves. The textured fingertips provide a secure grip, making these gloves suitable for tough jobs. At 12" long, they cover your wrist and forearm for added safety. Thickster® gloves are great for janitorial work, marine tasks, plumbing, painting, automotive work, and handling chemicals. Choose Thickster® gloves for confidence and comfort in demanding jobs!



### INDUSTRIES & APPLICATIONS

- Automotive
- Chemical
- Dental
- Electronics
- Food Industries
- Healthcare
- Janitorial
- Manufacturing
- Mining
- Oil and Gas
- Plumbing
- Recycling and Waste Disposal
- Transportation
- Utilities
- Warehousing
- Painting

### PRODUCT FEATURES

- Powder-free exam grade or powdered latex
- Thickness: 14 mil
- 12" length for added protection
- Beaded cuff for easy donning
- Fully textured for superior grip
- Strong yet soft and flexible
- Powder-free: Made from component materials which comply with FDA regulations for food contact. 21 CFR 170-199





## STANDARDS & PERFORMANCE

<b>Powder-Free FDA Regulations</b>	21 CFR 170-199 Food contact
<b>Inspection Level</b>	AQL 10.0 Level S3
<b>Standards</b>	ASTM D3578
<b>Storage</b>	Keep out of direct sunlight; store in a cool and dry place.

## TECHNICAL DATA

<b>Material</b>	Latex
<b>Thickness</b>	14 mil
<b>Length</b>	12"
<b>Grade</b>	Industrial & Exam
<b>Powder</b>	Powdered & Powder-Free
<b>Color</b>	Blue
<b>Texture</b>	Fully
<b>Sizes</b>	M - 2XL
<b>Usage</b>	Single Use
<b>Country of Origin</b>	Malaysia

## PHYSICAL PROPERTIES

<b>Length (mm / inches)</b>	305 / 12
<b>Palm Thickness (mm/mil)</b>	0.32 / 12.59
<b>Finger Thickness (mm/mil)</b>	0.36 / 14.17
<b>Tensile Strength</b>	22 MPa (min)
<b>Tensile Strength (Aged)</b>	18 MPa (min)
<b>Elongation</b>	750%
<b>Elongation (Aged)</b>	680%

## ITEM INFORMATION

### Powder-Free

#### Box of 50

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF
6602-20	Medium			2.1			22.18	781311660229	10781311660226
6603-20	Large	50 Gloves	10.5 x 5.25 x 3	2.41	10 Boxes	16 x 11.5 x 11.75	25.7	781311660328	10781311660325
6604-20	X-Large	per box		2.57	per case		27	781311660427	10781311660424
6605-20	2X-Large			2.67			28.13	781311660526	10781311660523

### Powdered

#### Box of 50

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF
6602	Medium			2.2			23.4	781311066021	10781311066028
6603	Large	50 Gloves	10.5 x 5.5 x 3.25	2.51	10 Boxes	16 x 11.38 x 11.5	26.5	781311066038	10781311066035
6604	X-Large	per box		2.6	per case		27.94	781311066045	10781311066042
6605	2X-Large			2.57			29.81	781311066052	10781311066059

### Pack of 4

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Inner Qty	Inner Dims (in)	Inner Wt (lb)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF	ITF
6602-22	S/M	4 Gloves	8 x 4 x .05	0.15	25 Packs	10 x 4 x 8	4.35	50 Packs	10 x 10 x 9	9.6	781311701021	30781311701022	50781311701026
6603-22	L/XL	per pack		0.15	per inner	10 x 4 x 7	5.1	per case		10.8	781311663220	10781311663227	20781311663224

Read and Understand Instructions and Warnings before using this product.

Raven Safety makes no recommendations as to the use of this product. All safety equipment should be selected in conjunction with an OSHA approved safety program that evaluates your specific hazards. Improper selection and/or use of safety products may cause serious injury.

