




# THICKSTER®



## THICKSTER EXAM GRADE LATEX DISPOSABLE GLOVES

Thickster® gloves offer excellent protection and performance. 14 mils thick, more than double the thickness of standard latex gloves. The textured fingertips provide a secure grip, making these gloves suitable for tough jobs. At 12" long, they cover your wrist and forearm for added safety. Thickster® gloves are great for janitorial work, marine tasks, plumbing, painting, automotive work, and handling chemicals. Choose Thickster® gloves for confidence and comfort in demanding jobs!



### INDUSTRIES & APPLICATIONS

- Automotive
- Chemical
- Dental
- Electronics
- Food Industries
- Healthcare
- Janitorial
- Manufacturing
- Mining
- Oil and Gas
- Plumbing
- Recycling and Waste Disposal
- Transportation
- Utilities
- Warehousing
- Painting

### PRODUCT FEATURES

- Powder-free exam grade or powdered latex
- Thickness: 14 mil
- 12" length for added protection
- Beaded cuff for easy donning
- Fully textured for superior grip
- Strong yet soft and flexible
- Powder-free: Made from component materials which comply with FDA regulations for food contact. 21 CFR 170-199





## THICKSTER™

### STANDARDS & PERFORMANCE

|  |  |
|--|--|
| <b>Powder-Free<br/>FDA Regulations</b> | 21 CFR 170-199 Food contact                                    |
| <b>Inspection Level</b>                | AQL 10.0 Level S3  |
| <b>Standards</b>                       | ASTM D3578   |
| <b>Storage</b>                         | Keep out of direct sunlight;<br>store in a cool and dry place. |

### TECHNICAL DATA

|                          |                        |
|--------------------------|------------------------|
| <b>Material</b>          | Latex                  |
| <b>Thickness</b>         | 14 mil                 |
| <b>Length</b>            | 12"                    |
| <b>Grade</b>             | Industrial & Exam      |
| <b>Powder</b>            | Powdered & Powder-Free |
| <b>Color</b>             | Blue                   |
| <b>Texture</b>           | Fully                  |
| <b>Sizes</b>             | M - 2XL                |
| <b>Usage</b>             | Single Use             |
| <b>Country of Origin</b> | Malaysia               |

### PHYSICAL PROPERTIES

|                                  |              |
|----------------------------------|--------------|
| <b>Length (mm / inches)</b>      | 305 / 12     |
| <b>Palm Thickness (mm/mil)</b>   | 0.32 / 12.59 |
| <b>Finger Thickness (mm/mil)</b> | 0.36 / 14.17 |
| <b>Tensile Strength</b>          | 22 MPa (min) |
| <b>Tensile Strength (Aged)</b>   | 18 MPa (min) |
| <b>Elongation</b>                | 750%         |
| <b>Elongation (Aged)</b>         | 680%         |

### ITEM INFORMATION

#### Powder-Free

##### Box of 50

| Item No. | Size     | Retail Qty           | Retail Dims (in) | Retail Wt (lbs) | Case Qty             | Case Dims (in)    | Case Wt (lbs) | UPC          | ITF            |
|----------|----------|----------------------|------------------|-----------------|----------------------|-------------------|---------------|--------------|----------------|
| 6602-20  | Medium   | 50 Gloves<br>per box | 10.5 x 5.25 x 3  | 2.1             | 10 Boxes<br>per case | 16 x 11.5 x 11.75 | 22.18         | 781311660229 | 10781311660226 |
| 6603-20  | Large    |                      |                  | 2.41            |                      |                   | 25.7          | 781311660328 | 10781311660325 |
| 6604-20  | X-Large  |                      |                  | 2.57            |                      |                   | 27            | 781311660427 | 10781311660424 |
| 6605-20  | 2X-Large |                      |                  | 2.67            |                      |                   | 28.13         | 781311660526 | 10781311660523 |

#### Powdered

##### Box of 50

| Item No. | Size     | Retail Qty           | Retail Dims (in)  | Retail Wt (lbs) | Case Qty             | Case Dims (in)    | Case Wt (lbs) | UPC          | ITF            |
|----------|----------|----------------------|-------------------|-----------------|----------------------|-------------------|---------------|--------------|----------------|
| 6602     | Medium   | 50 Gloves<br>per box | 10.5 x 5.5 x 3.25 | 2.2             | 10 Boxes<br>per case | 16 x 11.38 x 11.5 | 23.4          | 781311066021 | 10781311066028 |
| 6603     | Large    |                      |                   | 2.51            |                      |                   | 26.5          | 781311066038 | 10781311066035 |
| 6604     | X-Large  |                      |                   | 2.6             |                      |                   | 27.94         | 781311066045 | 10781311066042 |
| 6605     | 2X-Large |                      |                   | 2.57            |                      |                   | 29.81         | 781311066052 | 10781311066059 |

Read and Understand Instructions and Warnings before using this product.

Raven Safety makes no recommendations as to the use of this product. All safety equipment should be selected in conjunction with an OSHA approved safety program that evaluates your specific hazards. Improper selection and/or use of safety products may cause serious injury.

