



# **Industries and Applications**

- Automotive
  - Chemical
- Dental
- Electronics
- Food Industries
- Healthcare

- Janitorial
- Mining
- Oil and Gas
- Plumbing Recycling and Waste Disposal
- Transportation
- Utilities
- Warehousing Manufacturing Painting

### **Product Features**

- Powder-free or powdered exam grade latex
- Thickness: 14 mil
- 12" length for added protection
- Beaded cuff for easy donning
- Fully textured for superior grip
- Strong yet soft and flexible
- Powder-free: Made from component materials which comply with FDA regulations for food contact. 21 CFR 170-199

EXTENDED

# **THICKSTER**









# **Tech Specs**

Material Latex
Color Blue
Thickness 14 mil
Grade Exam
Length 11.81"
Country of Origin Malaysia
Available Sizes M - 2XL
Usage Single Use

# **Physical Properties**

Length (mm / inches) 300 / 11.81 Cuff Thickness (mm/mil) 0.2 / 7.87 Palm Thickness (mm/mil) 0.32 / 12.6 Finger Thickness (mm/mil) 0.36 / 14.17 Tensile Strength 22 MPa (min) Tensile Strength (Aged) 18 MPa (min) Elongation 750% **Elongation (Aged)** 680%

Width

Storage Keep out of direct sunlight;

store in a cool and dry place.

### **Item Information**

| Powder-Free<br>Box of 50<br>Item No.       | Size             | Retail Qty                 | Retail Dims (in)                    | Retail Wt (lbs)   | Case Qty                   | Case Dims (in)      | Case Wt (lbs)        | UPC               | ITF10          |
|--|------------------|----------------------------|-------------------------------------|-------------------|----------------------------|---------------------|----------------------|-------------------|----------------|
| 6602-20                                    | Medium           |                            |                                     | 2.1               |                            |                     | 22.18                | 781311660229      | 10781311660226 |
| 6603-20                                    | Large            | 50 Gloves<br>per box       | 10.5 x 5.25 x 3                     | 2.41              | 10 Boxes<br>per case       | 16 x 11.5 x 11.75   | 25.7                 | 781311660328      | 10781311660325 |
| 6604-20                                    | X-Large          |                            |                                     | 2.57              |                            |                     | 27                   | 781311660427      | 10781311660424 |
| 6605-20                                    | 2X-Large         |                            |                                     | 2.67              |                            |                     | 28.13                | 781311660526      | 10781311660523 |
| Powdered<br>Box of 50                      |                  |                            |                                     |                   |                            |                     |                      |                   |                |
| Item No.                                   | Size             | Retail Qty                 | Retail Dims (in)                    | Retail Wt (lbs)   | Case Qty                   | Case Dims (in)      | Case Wt (lbs)        | UPC               | ITF10          |
| 6602                                       | Medium           |                            |                                     | 2.2               |                            |                     | 23.4                 | 781311066021      | 10781311066028 |
| 6603                                       | Large            | 50 Gloves<br>per box       | 10.5 x 5.5 x 3.25                   | 2.515             | 12 Boxes<br>per case       | 16 x 11.38 x 11.5   | 24.3                 | 781311066038      | 10781311066035 |
| 6604                                       | X-Large          |                            |                                     | 2.6               |                            |                     | 27.935               | 781311066045      | 10781311066042 |
| 6605                                       | 2X-Large         |                            |                                     | 2.57              |                            |                     | 29.81                | 781311066052      | 10781311066059 |
| Powder-Free<br>Pack of 2 Pairs<br>Item No. | Size Retail Q    | ty Retail Dims (in)        | Retail<br>Wt (lbs) Inner Qty        | Inner Dims (in)   | Inner<br>Wt (lbs) Case Qty | Case Dims (in)      | Case<br>Wt (lbs) UPC | ITF10             | ITF20          |
| 6603-22                                    | Large 4 Gloves   | 4 GIOVES 9 5 x 5 88 x 0 38 | 0.205 100 Gloves<br>0.195 per inner | 15.25 x 12 x 3.75 | 5.21 4 Inners              | 21 13 8 8 25 8 11 5 | 21.925 7813116632    | 20 10781311663227 | 20781311663224 |
| 6604-22                                    | X-Large 4 Gloves |                            |                                     | 10.20 x 12 x 0.70 | per case                   |                     | 22.56 7813116642     | 27 10781311664224 | 20781311664221 |

These gloves are not intended for use as a chemical barrier for protection against all chemical compounds. Do not use if contaminated. Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for uses. SAS Safety assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

SAS Safety Corp. makes no recommendations as to the use of this product. All safety equipment should be selected in conjunction with an OSHA approved safety program that evaluates your specific hazards. Improper selection and/or use of safety products may cause serious injury.

Scan for more info!

