



THICKSTER



Thickster Latex Exam Grade Disposable Gloves

Thickster exam grade latex disposable gloves are 14 mil thick for extreme durability and dexterity. Made with a 12" length for added protection. Ideal for janitorial, marine, plumbing, painting, automotive and chemical resistance. Available in powdered or powder-free.



Industries and Applications

- Automotive
- Chemical
- Dental
- Electronics
- Food Industries
- Healthcare
- Janitorial
- Manufacturing
- Mining
- Oil and Gas
- Plumbing
- Recycling and Waste Disposal
- Transportation
- Utilities
- Warehousing
- Painting

Product Features

- Powder-free or powdered exam grade latex
- 12" length for added protection
- Beaded cuff for easy donning
- Fully textured for superior grip
- Single use only
- Strong yet soft and flexible
- Thickness: 14 mil



Tech Specs

Material	Latex
Color	Blue
Thickness	14 mil
Grade	Exam
Length	11.81"
Country of Origin	Malaysia
Available Sizes	M - 2XL
Usage	Single Use

Physical Properties

Length (mm / inches)	300 / 11.81
Cuff Thickness (mm/mil)	0.2 / 7.87
Palm Thickness (mm/mil)	0.32 / 12.6
Finger Thickness (mm/mil)	0.36 / 14.17
Tensile Strength	22 MPa (min)
Tensile Strength (Aged)	18 MPa (min)
Elongation	750%
Elongation (Aged)	680%

Width	
Small	97 ± 5 mm
Medium	107 ± 5 mm
Large	117 ± 5 mm
X-Large	115 ± 5 mm
2X-Large	120 ± 5 mm

Storage Keep out of direct sunlight; store in a cool and dry place.

Item Information

Powder-Free

Box of 50

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF10
6602-20	Medium			2.1			22.18	781311660229	10781311660226
6603-20	Large	50 Gloves	10.5 x 5.25 x 3	2.41	10 Boxes	16 x 11.5 x 11.75	25.7	781311660328	10781311660325
6604-20	X-Large	per box		2.57	per case		27	781311660427	10781311660424
6605-20	2X-Large			2.67			28.13	781311660526	10781311660523

Powdered

Box of 50

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF10
6602	Medium			2.2			23.4	781311066021	10781311066028
6603	Large	50 Gloves	10.5 x 5.5 x 3.25	2.515	12 Boxes	16 x 11.38 x 11.5	24.3	781311066038	10781311066035
6604	X-Large	per box		2.6	per case		27.935	781311066045	10781311066042
6605	2X-Large			2.57			29.81	781311066052	10781311066059

Powder-Free

Pack of 2 Pairs

Item No.	Size	Retail Qty	Retail Dims (in)	Retail Wt (lbs)	Inner Qty	Inner Dims (in)	Inner Wt (lbs)	Case Qty	Case Dims (in)	Case Wt (lbs)	UPC	ITF10	ITF20
6603-22	Large	4 Gloves	9.5 x 5.88 x 0.38	0.205	100 Gloves	15.25 x 12 x 3.75	5.21	4 Inners	21.13 x 8.25 x 11.5	21.925	781311663220	10781311663227	20781311663224
6604-22	X-Large			0.195	per inner			per case		22.56	781311664227	10781311664224	20781311664221

These gloves are not intended for use as a chemical barrier for protection against all chemical compounds. Do not use if contaminated. Components used in making these gloves may cause allergic reactions in some users. Follow your institution's policies for uses. SAS Safety assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

SAS Safety Corp. makes no recommendations as to the use of this product. All safety equipment should be selected in conjunction with an OSHA approved safety program that evaluates your specific hazards. Improper selection and/or use of safety products may cause serious injury.

Scan for more info!

